INTRAOCULAR LENS/DEVICE DESCRIPTION: The Focus Acrylic™ UV-Absorbing Posterior Chamber intraocular foldable lenses are precision made optical devices, designed for capsular bag implantation. The optical portion is biconvex with straight edges and made from high refractive index UV absorbing hydrophobic acrylate and has the capability of being folded. The Focus Acrylic Model 600 & 601 are transparent in color, Model 602 is yellow in color and it permanently blocks blue light while Model 603 (eclipse) is photochromic in nature and this property is activated only under the presence of UV light, under which the lens becomes yellow in color and blocks blue light.

Model	Optic Dia/Type	Overall Dia	Haptic	Diopter	Increments
600	6.0 mm Spheric Clear	12.5mm	Blue PVDF	0.0 to +9.0 D +10.0 to +30.0 D	1.0D 0.5 D
601	6.0mm Aspheric Clear	13.0mm	Acrylate	0.0 to +9.0 D +10.0 to +30.0 D +31.0 to +34.0 D	1.0D 0.5D 1.0D
602	6.0mm Aspheric Yellow	13.0mm	Acrylate	-7.0 to +9.0 D +10.0 to +30.0 D +31.0 to +34.0 D	1.0D 0.5D 1.0D
603	6.0mm Aspheric Photochromic	13.0mm	Acrylate	0.0 to +9.0 D +10.0 to +30.0 D +31.0 to +34.0 D	1.0 D 0.5 D 1.0D

INDICATIONS: Visual correction of aphakia in adults in whom a cataractous lens has been removed by phacoemulsification. The lens is intended for placement in the capsular bag.

PRECAUTIONS: Do not resterilize the lens by any method. Contact eyePx for return policy if lens becomes non-sterile. Use only sterile intraocular irrigation solutions and viscoelastics to rinse and/or lubricate the lens. Do not store the lens in direct sun light, or at temperatures over 50°C (122°F). Do not autoclave the lens. Reuse and/or re-sterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

WARNINGS: 1. Physicians considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio: a) Recurrent severe anterior or posterior segment inflammation or Uveitis. b) The use of Silicone oils in patients with current Vitreoretinal Disease or those at high risk for future disease that may require silicone oil as part of therapy should be reconsidered (Apple, et al. 1997) c) Surgical difficulties at the time of cataract extraction which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss. d) A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible. e) Circumstances that would result in damage to the endothelium during implantation. f) Suspected microbial infection. g) Children under the age of 2 years are not suitable candidates for intraocular lenses. h) Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.

COMPLICATIONS/ADVERSE EVENTS: The complications related to the implantation of any intraocular lens are essentially the same as for cataract surgery and may require secondary surgical intervention. Complications seen with the same type of IOL may include but are not limited to: Corneal Edema, Iritis, Hyphema, Macular Edema, Pupillary Block, Secondary Glaucoma, Cyclitic Membrane, Vitritis, Endophthalmitis, Retinal Detachment, Lens Dislocation. Adverse events seen with the same type of IOL may include but are not limited to: Hypopyon, Intraocular Infection, Acute Corneal Decompensation, Secondary Surgical Intervention: a) Lens replacement/removal, b) Retinal Detachment Repair, c) Repositioning of lens, d) Iridectomy, e) Vitrectomy, f) Wound Repair Leak, g) Photocoagulation, h) Removal of Residual Cortex material, i) Anterior Capsulotomy.

DETAILED DEVICE DESCRIPTION:

OPTICS

Material: Ultraviolet-absorbing high refractive index Hydrophobic Acrylate
Configuration: Biconvex Aspheric or Spheric with straight edge
Power: -7.0 through +34.0 diopter
UV cutoff at 10% T: 384 nm (+12.0 diopter lens) 387 nm (+27.0 diopter lens)

HAPTICS:	(Three Piece 600)	(One Piece 601, 602 & 603)
Material:	Polyvinylidene fluoride (PVDF)	Acrylate
Color:	Blue	Clear, Yellow or Photochromic
Configuration:	Modified-C	Modified-L
Angulation:	5 Degrees	Planar

DIRECTIONS FOR USE: 1. Examine the label on the lens package for proper lens model, power and expiration date. 2. Open the pouch and remove the lens in a sterile environment. 3. Examine the lens thoroughly to ensure that particles have not become attached to it and examine the lens optical surfaces for other defects. 4. The following folders and viscoelastics are recommended for lens implantation and handling. Folders: EL09 (Epsilon). Insertors: EL10 or EL18 (Epsilon). Viscoelastics include Healon®, Amvisc™, Duovisc™, Biolon™, Staar Visc™ II, Ophthalmic™ and/or Dispersa. 5. Folding may be performed at the typical orientations (6-12, 3-9, or 4-10) 6. Hydration is not required. 7. The use of viscoelastics during folding and injection is required. 8. The following injector cartridge systems have been qualified: Epsilon Injector Model EL-22 with Operaid cartridge Models OD502/OD522 by Ophtec. Medical LP604235C/LI604215 (single Piece) LC604220M/LI604210 (three piece). Disposable Injector OD655 (Three Piece lens) OD 665 (Single Piece Lens) by Ophtec. Disposable Medical Accuject Model ASB2200 (Single Piece).

Caution: Do not use lens if the package has been damaged. The sterility of the lens may have been compromised. Use of other non-validated cartridges may leave residual coating on the lens. Use of other non-validated folders/insertors/viscoelastics may damage the lens. When using folders the lenses that have been folded should be released into the eye within 5 minutes from the time of folding. When using injector the lens should be injected within two minutes after loading.

LENS POWER CALCULATIONS: The Physician should determine preoperatively the power of the lens to be implanted. Lens power calculation methods are described in the following references:

Dr. Haigis web site for the User Group for Laser Interference Biometry at the U. of Wuerzburg, <http://www.augenklinik.uni-wuerzburg.de/ulib/index.htm>.

Holladay, J.T., Musgrove, K.H., Prager, T.C., Lewis, J.W., Chandler, T.Y., and Ruiz, R.S., "A three-part system for refining intraocular lens power calculations," *J. Cataract Refract. Surg.*, Vol. 14, pp.17-24, 1988.

Holladay, J.T., "Standardizing Constants for Ultrasonic Biometry, Keratometry and Intraocular Lens Power Calculations"; *J. Cataract Refract. Surg.*, Vol.23, pp.1356-1370, 1997.

Physicians requiring additional information on lens power calculations may contact eyePx.

REPORTING: Adverse events and/or potentially sight-threatening complications that may be reasonably regarded as lens related and that were not previously expected in nature, severity or degree of incidence should be reported to eyePx. Irvine, CA USA. Fax 1-909-984-4671 or by e.mail: info@eyepx.com This information is being requested from all implant surgeons in order to document potential long term effects of intraocular lens implantation. Outside of the U.S contact your local company representative.

HOW SUPPLIED/ EXPIRATION DATE: Each IOL is supplied sterile, in dry form, in a lens container sealed within a single sterile pouch pack. The package is sterilized by Gamma or Hydrogen Peroxide. Open only under sterile conditions. Sterility is guaranteed unless the pouch is damaged or otherwise compromised. Expiration date and sterilization method is indicated on the outside of the lens box. The lens should not be implanted after the indicated sterility expiration date. Do not store the lens at temperatures below -20°C (-4°F) or above 50°C (122°F)

RETURN/EXCHANGE POLICY: Contact eyePx. for return lens policy. Return the lens with proper identification and reason for return. Label return as a biohazard. Lenses with past sterility expiration date should also be returned with prior authorization.

DISCLAIMER OF LIABILITY: eyePx accepts no liability for the choice of method or technique to implant the lens or for the choice of the lens for a particular patient or patient's condition.

REFERENCES: Apple, D.J., et al "Silicone Oil Adhesion to Intraocular Lenses: An experimental Study Comparing Various Biomaterials". *J.Cataract Refract Surg.* Vol 23, May 1997.

Manufacturer: eyePx LLC. 13 Avenue de Montrouge E. 6. R. Normandie 92340 Bourg la Reine, France

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Holladay, J.T., Musgrave, K.H., Prager, T.C., Lewis, J.W., Chandler, T.Y., and Ruiz, R.S., "A three-part system for refining intraocular lens power calculations," *J. Cataract Refract. Surg.*, Vol. 14, pp. 17-24, 1988.

Holladay, J.T., "Standardizing Constants for Ultrasonic Biometry, Keratometry and Intraocular Lens Power Calculations"; *J. Cataract Refract. Surg.*, Vol. 23, pp. 1356-1370, 1997.

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REFERENCES: Apple, D.J., et al "Silicone Oil Adhesion to Intraocular Lenses: An experimental Study Comparing Various Biomaterials". *J. Cataract Refract Surg.* Vol 23, Mai 1997.

Fabricant: eyePx LLC. 13 Avenue de Montrouge E. 6. R. Normandie 92340 Bourg la Reine, France

EC	REP
----	-----

 Dr. Arezi Beldjenna Phd. 13 Avenue de Montrouge E6. R. Normandie 92340 Bourg la Reine, France.

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Fabricante: eyePx LLC. 13 Avenue de Montrouge E. 6. R. Normandie 92340 Bourg la Reine, France